	•						
	1			_	-	-	
User Needs Intended Uses	Design Inputs (Technical Requirement)	Acceptar Ideal/Goal	nce Criteria Limit	Design Outputs	Essential Design Outputs (Check if essential)	Design Verification (Outputs = Inputs)	Design Validation (User Needs Met)
Customer Requirements (e.g. Functional)							
Pipetting volume range	10-500 uL	1-1000 uL	10-500 uL	Test usability of various sized pipettes (1-20 uL to 200-100 uL sized pipettes) to determine usability of multiple sizes Record masses of varying increments of the size of t			
Accurate pipetting in terms of liquid volume delivery	2 uL increment accuracy (tolerance)	1 uL increment	2 uL increment	liquids to determine the volume pipetted. Test with low and high viscosity liquids.			
volume delivery	(tolerance)	Less than	Less than 850x850x600	Measure LENGTH x			
Must fit on a lab bench User-friendly interface	800x800x550 mm Preset protocols for PCR and calibration	750x750x500 mm	mm	WIDTH x HEIGHT Determine effectiveness of protocalls by comparing to desired protocall output			Conduct survey for GUI satisfaction
Transportable by two individuals throughout a lab	Approximately 30 pounds. Not physcally awkward to carry. Transportable by 2 individuals.	Less than 25 pounds or transportable by a single person	Less than 40 pounds or transportable by two people Programmable by person	Weigh device. Determine how many individuals necessary to transport a short distance (between lab spaces)			
Programmable by intermediate in arduino code	Uses basic arduino code	Programmable by anyone in lab	with significant arduino code knowledge	Determine difficulty of code. Determine whether all			Have user write arduino code
Must fit multiple components	plates (common in lab), pipette tips, waste disposal	Compatible with all components in lab	Uses fewer than 3 external components	components stated could be used Use a nanodrop or other protein/nucleic acid			
Must avoid contaminating samples with external proteins, nucleic acids, or degratory enzymes	Comparable RNA, DNA, protein concentrations to assays run by hand			concentration detecting device and compare statistically the differente to manual pipetting			
Must save time/resources	Time spent accomplishing a task using the device is less than if not using the device (should take less than ~20 minutes per pcr plate) After running a programmed task	Maximum time per plate: 10 minutes	Maximum time per plate: 30 minutes	Measure time elapsed during single pipetting task (1 plate)			
Does not require supervision	on the device, the user does not need to supervise the device while it completes the task Repeated positioning error is	No individuals needed to supervise machine after selecting protocol	Individual needed for supervision for first and last 5 minutes of process				
Must be at least as accurate as by hand in terms of positioning	(percent deviation from human error)	on protocol	Makes less errors than a person Trained individual	quantify errors			
Manual for use	Documentation for protocls that lab novices can follow	Any individual can use manual to use basic protocols	needed to use basic protocols	Determine level of training necessary for use of manual			
Product Performance Requirements (e.g.N							
Must sustain repeated use (maintain accuracy) - calibration (self-calibration)	Lifespan of 10 years Calibration/accuracy verification	Lifespan of 50 years	Lifespan of 7 years				
Method to verify accuracy	mode that the user can run						
Powered by standard wall outlet Interfaces with Other Systems	120V power supply						
Potential to interface or include other lab equipment (centrifuge)							
Program must document actions, experiments	Count number of pipetted wells, document media/fluid used						
Compatible with current technology (pipettes, tips, etc.)	Existing tips and pipettes and plates used at sci tech						
Other							
Name Edmund Tamas Takata	Signature Edmund Tamas Takata	Date 10/11/2017					
Nickolas Hartman	Nickolas Hartman	10/11/2017					